



Product Information
Version 1.0

ZEISS ApoTome.2

Optical Sections in Fluorescence Imaging



We make it visible.

Simply Brilliant: Perfect Optical Section Thickness for All Magnifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › Technology and Details
- › Service

Create optical sections of your fluorescent samples – free of scattered light. With structured illumination, you know that only the focal plane appears in your image: ApoTome.2 recognizes the magnification and moves the appropriate grid into the beampath. The system then calculates your optical section from three images with different grid positions without time lag. It's a totally reliable way to prevent scattered out-of-focus light, even in your thicker specimens. Yet your system remains just as easy to operate as always. You get images with high contrast in the best possible resolution – simply brilliant optical sections.



Simpler. More Intelligent. More Integrated.

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Perfect Images – with All Magnifications



Because your applications need different objectives, you need a system that gives you the best resolution for each one. ApoTome.2 automatically uses the right grid for your objective, selecting from three grids with different frequencies. With a defined optical section thickness in the region of a Rayleigh unit, the image is simply brilliant.

Optimum Results – Free Choice of Light Source and Dyes

From conventional HBO illumination to adjustment-free metal-halide lamp HXP 120 C to Colibri.2, the LED illumination source that is gentle on your samples: with ApoTome.2 you use exactly the light you need. ApoTome.2 also gives you the choice of fluorophores. Whether you work with DAPI, FITC, Rhodamin, Cy5 or with vital dyes such as GFP or mRFP, it's your decision, not the technology's. Just change the filter and your system automatically moves the grid to the correct position. From DAPI to Cy5, you get perfect optical sections for multi-channel imaging.

Brilliant Images – Even with Thick Specimens

Your optical section thickness is close to one Rayleigh unit, a value that stands for high axial resolution with a good signal-to-noise ratio. ApoTome.2 increases the resolution in Z direction compared to conventional fluorescence microscopy: you obtain brilliant optical sections that allow 3D-rendering, even from thick specimens.



Rat, hippocampus, triple fluorescence, maximum-projection of 3D image-stack. Objective: Plan-APOCHROMAT 63x/1.4 E. Fuchs & S. Bauch, DPZ Göttingen, Germany


Three Grids Deliver Optimal Optical Section Thickness

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Figure A:

Acquisition with conventional epifluorescence illumination

Emission light from areas outside of the focal plane is detected. Contrast and resolution are reduced, depending on thickness of specimen.



C. elegans, whole mount, green: GFP, blue: DAPI
Objective: Plan-APOCHROMAT 20x/0.8
Prof. Schnabel, TU Braunschweig, Germany




Figures B – D:

Optical sections with different thickness

No matter which magnification you are using – ApoTome.2 automatically places the right grid in the beampath of your microscope.

Reduction of unwanted background fluorescence increases with the grid frequency and the optical sections become thinner.

Structures from outside of the focal plane are suppressed (Fig. B, C and D). This improves contrast and resolution of the optical section. “Low grid” delivers the optimal section thickness in our example (Fig. D). Images of this type are particularly suitable for 3D analyses and the processing of your image data with rendering software.



Your Insight into the Technology Behind It


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

ZEISS ApoTome.2

Brings You Structured Illumination

ApoTome.2 projects a grid structure into the focal plane of your specimen, then moves it into three positions using a scanning mechanism.

At each grid position, ApoTome.2 automatically acquires a digital image. The system processes the three images into one optical section with improved contrast and increased resolution using a patented algorithm. The image that emerges is free from grid structures.



Schematic illustration of the grid projection.
A – C: raw images with different positions of grid
D: optical section through sample



Animation from www.zeiss.com/campus, © Mike Davidson, FSU, Tallahassee

ZEISS ApoTome.2 Grid in the Beampath

Fluorescence excitation light passes through two glass plates in the ApoTome.2 slider. When a grid structure is applied to the first glass plate, the grid pattern is “imprinted” in the excitation light. A scanning mechanism tilts the second glass plate and the image of the grid is laterally shifted in the focal plane of the specimen.

Tailored Precisely to Your Applications

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

ApoTome.2 is the cost-effective solution for creating optical sections with high contrast. Use this to your advantage in a wide range of applications from cell culture preparations via tissue sections to whole embryos.

| Typical Applications / Typical Specimens | Task | ZEISS ApoTome.2 offers |
|--|---|---|
| Cell Culture | 2D imaging | 2D single images possible |
| | Fast imaging of a 2D image | Optical section available online on the monitor |
| | Reliable detection of the marker even with strong background fluorescence | Automatic grid selection ensures optimum contrast with each objective |
| | Combination of multiple contrast techniques | Any combination of fluorescence channels, brightfield, DIC and phase contrast. Each fluorescence channel can be individually configured as an optical section or widefield image |
| Live Cell Imaging | Reduction of phototoxicity | Particularly low in combination with LED illumination and EMCCD cameras |
| | Time-lapse images | Depending on the exposure time, up to three images per second. Doubling of the frame rate with "burst mode" |
| Vibratome Sections, Histological Samples | 3D imaging | Automatic selection of the optimum grid for each objective |
| | Modification of the optical section thickness | Grid freely selectable depending on the specimen |
| | Penetration depth | Depending on the optical density of the tissue |
| | 3D reconstruction | Rendering of the image stack via integrated software function. Automatic transfer of the parameters of the individual fluorescence channels |
| | Quantitative analysis | Automatic calibration of the system: reproducible size measurements |
| Whole Mounts | 3D imaging | Multi Channel, Z Stack and Time Lapse, Deconvolution, images in raw data mode, 3D Rendering |
| | Large image areas | Automatic acquisition of large sections using Tiles & Positions |

ZEISS ApoTome.2 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service




Figure A: Conventional fluorescence





Figure B: Optical section

Drosophila neurons, blue: DAPI, green: GFP. Objective: Plan-APOCHROMAT 20 x/0.8.
Marta Koch, Molecular and Developmental Genetics, University of Leuven, Belgium



▶ [Click here to view this video on YouTube](#)

Figure C: *Drosophila* embryo, green: HRP, red: glia marker,
100 μ m Z-stack
C. Klämbt, Institute for Neurobiology, University of Münster,
Germany




▶ [Click here to view this video on YouTube](#)

Figure D: Mouse embryo, tissue section, green: GFP, red: Cy3
Objective: Plan APOCHROMAT 40 x/1.3 Oil
N. Büttner, T. Vogel, Centre for Anatomy, University of Göttingen, Germany

Your Flexible Choice of Components

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



1 Microscope

- Axio Observer, Axiovert 200 (inverted research microscope)
- Axio Imager.2, Axio Imager.Z1, Axio Imager.D1 (upright research microscope)
- Axio Zoom.V16 (Zoom microscope)
- Simple upgrading of existing systems

2 Objectives

Recommended objective classes with the highest level of image quality:

- C-APOCHROMAT
- Plan-APOCHROMAT
- EC Plan-NEOFLUAR

3 Illumination

- Colibri (LED)
- HXP 120 C (metal halide)
- HBO (mercury vapor lamp)
- XBO (xenon)

4 Cameras

- Recommended cameras with high dynamic range (thick samples: at least 1 : 2000; thin samples at least 1 : 1000; digitalization at least 12 bit)
- AxioCam MRm

- Alternatively, you can control these cameras: Photometrics CoolSnap HQ, Hamamatsu Orca ER2 (cameras with pixel size providing a sampling rate of < 5. This allows a pixel size of approx. 6.5 μm (Imaging in Neuroscience and Development 2005; Chapter 101; pp. 805–813)

5 Software

Recommended ZEN modules:

- Multi Channel, Z Stack, Time Lapse (imaging)
- Tiles & Positions (imaging with scanning table)
- 3D VisArt (rendering multidimensional image stacks)
- Image analysis modules such as Image Analysis, Colocalization



Expand Your Possibilities

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service

Deconvolution

Improve the image stacks you create with ApoTome.2 even more with deconvolution, using a patented algorithm for structured illumination:

- Acquire image stacks in raw data format – individual images are saved for the different grid positions.
- Switch between conventional fluorescence and optical section after image acquisition.
- Deconvolution processes the raw data with a special algorithm for structured illumination.
- Enjoy improved image quality, contrast, axial and lateral resolution.
- The efficient suppression of any existing noise improves recognition of the object structures.



Example image of yeast cells: (above) optical section, (below) result of deconvolution.

Literature:


L. H. Schaefer, D. Schuster & J. Schaffer, "Structured illumination microscopy: Artefact analysis and reduction utilizing a parameter optimization approach", *Journal of Microscopy*, Vol. 216, Pt 2 November 2004, pp. 165–174.

System Overview

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service

ApoTome.2 is Compatible with these Stands from ZEISS:

- Axioplan 2 imaging
(serial numbers: from 35 11 000001;
from 35 10 000001; from 35 02 000001)
- Axio Imager.D1 and
Axio Imager.Z1, Axio Imager.A2
- Axio Imager.M2
- Axio Imager.D2 and Axio Imager.Z2
- Axiovert 200M, Axio Observer.A1
- Axio Observer.D1 and Axio Observer.Z1
- Axio Zoom.V16



Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Grid Table: ApoTome.2 generates optical sections of a defined thickness (in Rayleigh units, RU and microns, μm) depending on wavelength, microscope and objective used.

| Data for the Use of Upright Microscopes, e.g. ZEISS Axio Imager | | | | | | | | |
|---|------|------|------------------------|--|-------------|----------|----------------|----------------|
| Objectives for Axio Imager | V | NA | Immersion | Grid/Section thickness @490nm [RU/ μm] | | | DAPI with FS34 | DAPI with FS49 |
| | | | | High grid | Medium grid | Low grid | | |
| EC Plan-NEOFLUAR | 10x | 0.3 | Air | 2.9/31.9 | 1.7/18.2 | 0.9/9.9 | Yes | Yes |
| EC Plan-NEOFLUAR | 20x | 0.5 | Air | 2.4/9.2 | 1.4/5.3 | 0.7/2.9 | Yes | Yes |
| EC Plan-NEOFLUAR | 40x | 0.75 | Air | 1.6/2.8 | 0.9/1.6 | 0.5/0.9 | Yes | Yes |
| EC Plan-NEOFLUAR | 40x | 1.3 | Oil | 2.5/2.2 | 1.4/1.2 | 0.8/0.7 | Yes | Yes |
| EC Plan-NEOFLUAR | 63x | 0.95 | Air | 1.0/1.1 | 0.6/0.7 | 0.4/0.4 | Yes | No |
| EC Plan-NEOFLUAR | 63x | 1.25 | Oil | 1.6/1.5 | 0.9/0.9 | 0.5/0.5 | Yes | Yes |
| EC Plan-NEOFLUAR | 100x | 1.3 | Oil | 1.0/0.9 | 0.6/0.5 | 0.4/0.3 | Yes | Yes |
| LCI Plan-NEOFLUAR | 25x | 0.8 | Oil, water or glycerin | 2.9/6.6 | 1.7/3.7 | 0.9/2.0 | Yes | Yes |
| LCI Plan-NEOFLUAR | 63x | 1.3 | Water or glycerin | 1.5/1.3 | 0.9/0.7 | 0.5/0.4 | Yes | Yes |
| Plan-APOCHROMAT | 10x | 0.45 | Air | 4.2/20.4 | 2.4/11.5 | 1.3/6.2 | Yes | Yes |
| Plan-APOCHROMAT | 20x | 0.8 | Air | 3.2/4.9 | 1.8/2.8 | 1.0/1.5 | Yes | Yes |
| Plan-APOCHROMAT | 40x | 0.95 | Air | 1.6/1.7 | 0.9/1.0 | 0.5/0.5 | Yes | Yes |
| Plan-APOCHROMAT | 40x | 1.3 | Oil | 2.5/2.2 | 1.4/1.2 | 0.8/0.7 | Yes | Yes |
| Plan-APOCHROMAT | 40x | 1.4 | Oil | 2.4/1.8 | 1.4/1.0 | 0.7/0.6 | Yes | Yes |
| Plan-APOCHROMAT | 63x | 1.4 | Oil | 1.6/1.2 | 0.9/0.7 | 0.5/0.4 | Yes | Yes |
| Plan-APOCHROMAT | 100x | 1.4 | Oil | 1.0/0.8 | 0.6/0.5 | 0.4/0.3 | Yes | Yes |
| LD LCI Plan-APOCHROMAT | 25x | 0.8 | Oil, water or glycerin | 2.9/6.6 | 1.7/3.7 | 0.9/2.0 | Yes | Yes |
| C-APOCHROMAT | 10x | 0.45 | Water | 4.2/20.4 | 2.4/11.5 | 1.3/6.2 | Yes | Yes |
| C-APOCHROMAT | 40x | 1.2 | Water | 2.2/2.0 | 1.2/1.1 | 0.7/0.6 | Yes | Yes |
| C-APOCHROMAT | 63x | 1.2 | Water | 1.4/1.3 | 0.8/0.7 | 0.5/0.4 | Yes | Yes |
| LD C-APOCHROMAT | 40x | 1.1 | Water | 2.2/2.3 | 1.2/1.3 | 0.7/0.7 | Yes | Yes |
| Plan-APOCHROMAT | 63x | 1.46 | Oil | 1.5/1.0 | 0.9/0.6 | 0.5/0.3 | Yes | Yes |
| Plan-FLUAR | 100x | 1.45 | Oil | 1.0/0.7 | 0.6/0.4 | 0.3/0.2 | No | No |
| Plan-APOCHROMAT | 100x | 1.46 | Oil | 1.0/0.7 | 0.6/0.4 | 0.3/0.2 | Yes | No |

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Data for the Use of Inverted Microscopes, e.g. ZEISS Axio Observer

| Objectives for Axio Imager | V | NA | Immersion | Grid/Section thickness @490nm [RU/μm] | | | DAPI with FS34 | DAPI with FS49 |
|----------------------------|------|------|------------------------|---------------------------------------|-------------|----------|----------------|----------------|
| | | | | High grid | Medium grid | Low grid | | |
| EC Plan-NEOFLUAR | 10x | 0.3 | Air | 2.9/31.5 | 1.7/18.5 | 0.9/9.8 | Yes | Yes |
| EC Plan-NEOFLUAR | 20x | 0.5 | Air | 2.3/9.0 | 1.4/5.4 | 0.7/2.9 | Yes | Yes |
| EC Plan-NEOFLUAR | 40x | 0.75 | Air | 1.6/2.7 | 0.9/1.6 | 0.5/0.9 | Yes | No |
| EC Plan-NEOFLUAR | 40x | 1.3 | Oil | 2.4/2.1 | 1.4/1.3 | 0.8/0.7 | Yes | Yes |
| EC Plan-NEOFLUAR | 63x | 0.95 | Air | 1.0/1.1 | 0.6/0.7 | 0.4/0.4 | Yes | Yes |
| EC Plan-NEOFLUAR | 63x | 1.25 | Oil | 1.6/1.5 | 0.9/0.9 | 0.5/0.5 | Yes | No |
| EC Plan-NEOFLUAR | 100x | 1.3 | Oil | 1.0/0.9 | 0.6/0.6 | 0.4/0.3 | Yes | No |
| LCI Plan-NEOFLUAR | 25x | 0.8 | Oil, water or glycerin | 2.9/6.5 | 1.7/3.8 | 0.9/2.0 | Yes | Yes |
| LCI Plan-NEOFLUAR | 63x | 1.3 | Water or glycerin | 1.5/1.3 | 0.9/0.8 | 0.5/0.4 | No | No |
| Plan-APOCHROMAT | 10x | 0.45 | Air | 4.2/20.2 | 2.4/11.7 | 1.3/6.1 | Yes | Yes |
| Plan-APOCHROMAT | 20x | 0.8 | Air | 3.1/4.8 | 1.8/2.8 | 1.0/1.5 | Yes | Yes |
| Plan-APOCHROMAT | 40x | 0.95 | Air | 1.6/1.7 | 0.9/1.0 | 0.5/0.5 | Yes | Yes |
| Plan-APOCHROMAT | 40x | 1.3 | Oil | 2.4/2.2 | 1.4/1.3 | 0.8/0.7 | Yes | Yes |
| Plan-APOCHROMAT | 40x | 1.4 | Oil | 2.4/1.8 | 1.4/1.1 | 0.7/0.6 | Yes | Yes |
| Plan-APOCHROMAT | 63x | 1.4 | Oil | 1.5/1.2 | 0.9/0.7 | 0.5/0.4 | Yes | Yes |
| Plan-APOCHROMAT | 100x | 1.4 | Oil | 1.0/0.8 | 0.6/0.5 | 0.4/0.3 | Yes | No |
| LD LCI Plan-APOCHROMAT | 25x | 0.8 | Oil, water or glycerin | 2.9/6.5 | 1.7/3.8 | 0.9/2.0 | Yes | Yes |
| C-APOCHROMAT | 10x | 0.45 | Water | 4.2/20.2 | 2.4/11.7 | 1.3/6.1 | Yes | Yes |
| C-APOCHROMAT | 40x | 1.2 | Water | 2.1/1.9 | 1.3/1.1 | 0.7/0.6 | Yes | Yes |
| C-APOCHROMAT | 63x | 1.2 | Water | 1.4/1.3 | 0.8/0.7 | 0.5/0.4 | Yes | Yes |
| LD C-APOCHROMAT | 40x | 1.1 | Water | 2.1/2.3 | 1.3/1.4 | 0.7/0.7 | Yes | Yes |
| Plan-APOCHROMAT | 63x | 1.46 | Oil | 1.5/1.0 | 0.9/0.6 | 0.5/0.3 | Yes | Yes |
| Plan-FLUAR | 100x | 1.45 | Oil | 1.0/0.7 | 0.6/0.4 | 0.3/0.2 | No | No |
| Plan-APOCHROMAT | 100x | 1.46 | Oil | 1.0/0.7 | 0.6/0.4 | 0.3/0.2 | Yes | No |

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

| Dimensions | (Width x Depth x Height) |
|--|--|
| ApoTome.2 slider for Axio Imager | Approx. 278 mm x 90 mm x 76 mm |
| ApoTome.2 slider for Axio Observer/Axiovert 200 | Approx. 295 mm x 90 mm x 78 mm |
| Control box ApoTome.2 | Approx. 255 mm x 220 mm x 96 mm |
| Operating Data | |
| Protection Class, Protection Type | I, IP 20 |
| Electrical Safety | According to DIN EN 61010-1 (IEC 61010-1) taking account of CSA and UL regulations |
| Overvoltage Category | II |
| Interference Suppression | In accordance with EN 55011 class B |
| Interference Resistance | In accordance with DIN EN 61326-1 |
| Supply Voltage | 100 to 240 V ±10%. No Adjustment of the supply voltage is required |
| Supply Frequency | 50 to 60 Hz |
| Power Consumption ApoTome.2 | Max. 50 VA |
| Fuses in Accordance with IEC 127 | |
| Control box ApoTome.2 | 2 A delayed-action/H/250 V, 5 x 20 mm |
| Grid Frequencies | |
| Axio Imager slider (transmission grid high/medium/low) | 5/9/17,5 lp/mm |
| Axio Observer/Axiovert 200 slider (transmission grid high/medium/low) | 10/17.5/35 lp/mm |
| Installation Conditions | |
| <p>The grid projection method used for the ApoTome.2 is sensitive to vibration, which can have various causes (including strong draughts). Vibrations are visible as streak artefacts in the resulting image. The microscope must therefore be set up so that it is exposed to as little vibration as possible on a vibration-damped table or on a suitable microscope base.</p> | |



Count on Service in the True Sense of the Word

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › Technology and Details
- › **Service**

Because the ZEISS microscope system is one of your most important tools, we make sure it is always ready to perform. What's more, we'll see to it that you are employing all the options that get the best from your microscope. You can choose from a range of service products, each delivered by highly qualified ZEISS specialists who will support you long beyond the purchase of your system. Our aim is to enable you to experience those special moments that inspire your work.

Repair. Maintain. Optimize.

Attain maximum uptime with your microscope. A ZEISS Protect Service Agreement lets you budget for operating costs, all the while reducing costly downtime and achieving the best results through the improved performance of your system. Choose from service agreements designed to give you a range of options and control levels. We'll work with you to select the service program that addresses your system needs and usage requirements, in line with your organization's standard practices.

Our service on-demand also brings you distinct advantages. ZEISS service staff will analyze issues at hand and resolve them – whether using remote maintenance software or working on site.

Enhance Your Microscope System.

Your ZEISS microscope system is designed for a variety of updates: open interfaces allow you to maintain a high technological level at all times. As a result you'll work more efficiently now, while extending the productive lifetime of your microscope as new update possibilities come on stream.



Profit from the optimized performance of your microscope system with services from ZEISS – now and for years to come.

>> www.zeiss.com/microservice

The moment your data change scientific minds.
This is the moment we work for.

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › Technology and Details
- › Service



// RECOGNITION
MADE BY ZEISS



Carl Zeiss Microscopy GmbH
07745 Jena, Germany
BioSciences
microscopy@zeiss.com
www.zeiss.com/apotome

**Instruments**
info@ctk-instruments.com
www.ctk-instruments.com
866-464-1005



We make it visible.